

Matrix ORientation and Texture EXplorer (MORTEX)

Instruction Manual

Overview

MORTEX is an ImageJ macro for automated processing and feature extraction of digital extracellular matrix (ECM) microscopy images obtained with fluorescence imaging of Picrosirius red-stained tissue sections.

For the latest version, please visit the repository at www.github.com/cjravensbergen/MORTEX.

MORTEX is available for public use under the terms of the *GNU General Public License (GPL) v3*.


*** Please follow best practice guidelines for image analysis and maintain fixed acquisition and processing parameters during experiments. ***

Interface

MORTEX offers a minimalistic user interface for running its analysis module. Upon starting the **MORTEX** ImageJ macro, users are prompted to input the image file extension.

Press `OK` and a new window will pop up, prompting the user to select an image input directory. A second window will pop up prompting the use to select a directory for the output results file (.txt).


MORTEX provides quantitative output for ten matrix and fiber parameters (see description table below). The macro is optimized for batch processing of multiple images, returning a dataframe with measurements to the user-specified output directory.



The next sections describe the recommended experimental setup and pre-processing steps of **MORTEX** and provides a description table of the ECM parameters measured.

Experimental Setup

MORTEX is designed to analyze microscopy images obtained from fluorescent imaging of Picrosirius red-stained tissue sections. The experimental setup required to generate input images for **MORTEX** is shown below. Fluorescent acquisition using FITC (for autofluorescence) and TRITC (for ECM visualization) filters is recommended. The use of regions of interest (ROI) is advised to reduce processing times.



MORTEX Pre-Processing Steps

MORTEX applies basic image pre-processing before image analysis. The sequential pre-processing steps are described below.

Split Channels

Multichannel images are split into single channel images, where the ECM fiber image (TRITC signal) is retained for downstream analysis.

8-bit conversion

Thresholding

Fiber images are thresholded using ImageJ's 'Auto Threshold' function with default method, black and white are ignored.

Denoising

Thresholded fiber images are denoised using ImageJ's 'Despeckle' function.

Skeletonization













Thresholded fiber images are skeletonized using ImageJ's 'Skeletonize' function.

ECM Parameter Description Table

ECM Parameter	Description
Tortuosity	Measures the complexity or winding nature of fiber paths within the ECM. Higher tortuosity indicates more twisting and turning paths.
Compactness	The ratio of the area of a fiber network to its perimeter, indicating the compactness or spread of the structure. A high ratio suggest a compact ECM.
Uniformity	Measures the texture of the ECM by analyzing similarity in pixel intensities of neighboring fibers. High uniformity indicates a uniform matrix with repetitive elements.
Intersection density	Quantifies the number of intersections or cross-points between fibers within a defined area, indicating the degree of connectivity and complexity of the network.
Fractal dimension	A measure of the ECM's complexity and how detail changes with scale. Higher values indicate more intricate, space-filling patterns that are self-similar across scales.
Fiber bundle density	Represents the concentration of fiber bundles in a given area, indicating the thickness and robustness of the ECM network. Higher density suggests more tightly packed fibers.
Dominant direction	Indicates the primary orientation or alignment of fibers within the ECM, represented as an angle. It shows the main structural direction of the ECM.
Branching density	Measures the number of branching points within the ECM, indicating the extent of network branching. High branching density suggests a complex, tree-like network structure.
Average fiber length	The mean length of individual fibers in the ECM. Longer fibers indicate elongated structures, while shorter fibers suggest more fragmented or crosslinked configurations.
Anisotropy index	Measures the degree of directional dependence in the ECM structure, showing how much the fibers are oriented in a particular direction. Higher anisotropy indicates more alignment.

Low value

High value



Contact & Contributions

Written by *Cor Ravensbergen* at the *Wilma Mesker lab*, Department of Surgery, Leiden University Medical Center, the Netherlands.

In collaboration with the *Light & Electron Microscopy Facility*, Cell & Chemical Biology, Leiden University Medical Center, the Netherlands.

Contact

c.j.ravensbergen@lumc.nl

www.github.com/cjravensbergen/MORTEX

Citation

Ravensbergen CJ, et al. ...

[\[PUBMED link\]](#)

Dependencies

- *Bio-Formats Importer*
- *GLCM2, version 1.0.1* (<https://github.com/miura/GLCM2>)
- *OrientationJ, version 2.0.5* (<https://bigwww.epfl.ch/demo/orientation/>)
- *AnalyzeSkeleton*